

Donévit Prostram

Micro magnetic and vibrant massage with heat against enlarged prostate trouble!

The device stimulates the prostate to recreate its more natural size and function. This by means of magnetism, vibration and heat to reduce the occurrence of fluid that prevents blood circulation in the prostate itself.



Enlarged Prostate problem and diseases bothers whole your life! From sleepless nights to huge stress.

Donévit Prostram

If you suffer from prostate problems, enlarged or sluggish prostate, then Donévit Prostram may be the best first aid against your discomfort.

Prostate problems can be remedied with Donévit Prostram developed in collaboration with the healthcare to "cure" and counteract enlarged Prostate diseases and problems.

Enlarged prostate can lead to serious illnesses and infections. The unique ability of the device that NOT adversely affect other bodily functions; have made it to be part of treatments in several hospitals, clinics and caregivers there this unit have got a given place as part in healthcare procedure against enlarged prostate.

The unit pulsates with variable frequency and low-effect magnetic field with heat (app 40°C) to ensure that the prostate regains more of the normal functions and its more natural size.

Recommended that it be used in parallel with the usual medical care. It is even easy to use in your own home in in the aftercare.

Donévit Prostram is intended for the treatment of patients with chronic prostatitis, internal effects on the prostate through three different methods, Hypertherm so called heat exposure, pulsing vibration and magnetism to get blood circulation started and help empty the prostate on "redundant" fluid formation and last but not least an anti-inflammatory function for to moderate Inflammations that can has due to catheters that are often offered while waiting for surgery.

Donévit Prostram will also improve your sex life. stronger self-esteem and better erection ability.


Benefits of Donévit Prostram prostate treatments:

- **reduced need of pee,**
- **less or no need of catheter**
- **minor infection risk**
- **better night's sleep**
- **less stress**
- **higher work effort**
- **and so on ...**



Contact:
sales@donevit.com
www.donevit.com
Mobil: +4673 59 59 200

Donévit Prostatam can in most cases help with enlarged prostate and prevent any inconvenience this may cause! Even after some treatments you can experience the first signs of a more comfortable life.




If you suffer from prostate, enlarged or sluggish prostate, Donévit Prostatam may be the best first aid to your problems.

Helps with the following:

- Chronic prostatitis
- Vascular disease (blood vessel disease and inflammation of the vessel wall)
- Dysfunction with erection (impotence)
- Prostate hyperplasia (benign)
- Dysfunction in microcirculation
- Inflammations (due to catheter handling)
- Strengthens the pelvic floor muscles
- Sluggish prostate (activates its functionality)
- Reduces irritation in the lower urinary tract
- Strengthens the muscle for sexual ejaculation at intercourse



Donévit Prostatam will be your help even if the problems are repeated. The device's efficacy increases several times when used in combination with conventional medical care. Donévit Prostatam does not replace medical care but is an excellent complement and a preventive tool. Hygiene protective rubber (condom) is not included without buying separately. You insert the device into the rectum and allow it to treat the prostate during approximately 25 minutes. 2 treatments a day initially then 1 time every other day until the problems decrease or decrease in scale. The device has a built-in timer that delivers signal when processing is complete.



Donévit Prostatam is suitable for professionals, therapists and home use!



Vi stödjer cancerfonden

Always read the manual before using the unit. Check out our YouTube Channel for info and news. The address can be found on our website: www.donevit.com

Donévit

Contact:
 sales@donevit.com
 www.donevit.com
 Mobil: +4673 59 59 200